

Parque Tecnológico de León. Parcela M-10.4  
24009 Armunia (León) SPAIN  
Telephone: +34 987 849 200

León 13/12/2018



Product Code:

AD0034

Name:

Napyradiomycin A; Napyradiomycin A1; WS-9558-A

## 1. General information

Molecular Formula:	C <sub>25</sub> H <sub>30</sub> Cl <sub>2</sub> O <sub>5</sub>
Molecular Weight:	480
CAS #:	103106-24-7
Long-term Storage:	Keep in dark, dry and cool place -20 °C
Assay (HPLC-UV, %)	83,79
Appearance	orange oily solid
Date of Last Analysis	09/12/2015
Retest date	08/12/2020

## 2. Activity according to literature

antibiotic, g-pos., estrogen receptor antagonist

## 3. Format

This compound is available in every format from 1 up to 5 mg. In case of needing more than 5 vials handling costs could be applied. Please feel free to send a request for amounts above 5 mg.

For some of our products grams might be immediately available

## 3. Advice and remarks

Whenever the container is opened for removal of aliquot portions of the substance, the person handling the substance must assure that the integrity of the substance is maintained and proper records of all its handlings are kept. Special care has to be taken to avoid any contamination or adulteration of the substance.

The pharmacological and toxicological properties of this product have not been fully investigated.

Analytical HPLC conditions by request.

Use caution when handling.

**For research use only. Not for human or drug use.**

Contact person: Laura Rodríguez Rodríguez

Phone: (+34) 987 849200 Email: l.rodriguez@biomar.co Web site: www.biomarmt.com