

1- Identification of the substance/preparation and of the company/undertaking

- **PRODUCT NAME:** Malformin A, from microbial source
- **PRODUCT NUMBER:** CH0057
- **CAS #:** 3022-92-2
- **SYNONYMS:** Malformin A1
- **MOLECULAR FORMULA:** C₂₃H₃₉N₅O₅S₂
- **MOLECULAR WEIGHT:** 529,7 g/mol
- **DETAILS OF THE MANUFACTURER:**

BIOMAR MICROBIAL TECHNOLOGIES
PARQUE TECNOLÓGICO DE LEÓN PARCELA M10.4
24009 ARMUNIA (LEÓN) SPAIN
TELF: +34987849200 FAX: +34987849203
e-mail: info@biomarmt.com
www.biomarmt.com

2- Composition/information on ingredients

- **Chemical characterization:** Substance
 - **CAS No. Description:** 3022-92-2, Malformin A
 - **Identification number(s):** Not applicable

3- Hazards identification

- **PHYSICAL/CHEMICAL HAZARDS:** Not available
- **HUMAN HEALTH HAZARDS:** Acute toxicity, Oral (Level 3)
Acute toxicity, Inhalation (Level 3)
Acute toxicity, Skin contact (Level 3)

4- First-aid measures

Description of first aid measures

- **General information:** Consult a physician. Show this safety data sheet to the doctor in attendance.

- **AFTER INHALATION:** If breathed in, move person into fresh air. If not breathing, give artificial respiration. Consult a physician.
- **AFTER INGESTION:** Never give anything by mouth to an unconscious person. Rinse mouth with water. Consult a physician.
- **AFTER SKIN CONTACT:** Wash off plenty of water while removing contaminated clothing and shoes. Consult a physician.
- **AFTER EYE CONTACT:** Rinse opened eye for several minutes under running water. Check for and remove contact lenses. If symptoms persist, consult a physician.
- **AFTER AGGRAVATING CONDITIONS:**

Indication of any immediate medical attention and special treatment needed.

No further relevant information available

5- Fire-fighting measures

- **FLAMMABILITY OF THE PRODUCT:** Not available
- **EXTINGUISHING MEDIA:** Water, dry chemical powder, foam, CO₂. Suitable extinguishing media for the surrounding fire should be used.
- **HAZARDOUS THERMAL (de) COMPOSITION PRODUCTS:** During a fire, irritating and highly toxic gases may be generated by thermal decomposition or combustion.
- **SPECIAL FIRE-FIGHTING PROCEDURES:** Wear self-contained breathing apparatus and protective clothing to prevent contact with skin and eyes.
- **PROTECTION OF FIRE-FIGHTERS:** Wear self-contained respiratory protective device. Do not inhale explosion gases or combustion gases. Wear fully protective suit.

6- Accidental release measures

- **PERSONAL PRECAUTIONS PROTECTIVE EQUIPMENT AND EMERGENCY PROCEDURES:**
Wear appropriate respirator, rubber boot and heavy rubber gloves. Scoop up and place in an appropriate container. Ventilate area and wash spill site after pickup is complete. Wash skin immediately with plenty of water.
- **ENVIRONMENTAL PRECAUTIONS:** No special measures required.
- **METHODS AND MATERIAL FOR CONTAINMENT AND CLEANING UP:** Absorb on sand or vermiculite and place in closed containers for disposal.

7- Handling and storage

- **HANDLING PRECAUTIONS FOR SAFE HANDLING:** Use caution when handling. Exposure to any chemical should be limited. Avoid prolonged exposure. Avoid all direct contact with product. Wear protective safety goggles. Wear chemical-resistant gloves. Wear protective clothing and boots. Ensure ventilation during use.

- **STORAGE CONDITIONS FOR SAFE STORAGE, INCLUDING ANY INCOMPATIBILITIES:** Store in cool place. Keep container tightly closed in a dry and well-ventilated place.

8- Exposures controls / personal protection

- **ADDITIONAL INFORMATION ABOUT DESIGN OF TECHNICAL FACILITIES:** No further data.
- **CONTROL PARAMETERS INGREDIENTS WITH LIMIT VALUES THAT REQUIRE MONITORING AT THE WORKPLACE:** Not required.
- **EXPOSURE CONTROLS:**
- **PERSONAL PROTECTIVE EQUIPMENT:**
 - **Skin and body:** Protective gloves. The glove material has to be impermeable and resistant to the product.
Due to missing tests no recommendation to the glove material can be given for the product. Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation
 - **Eyes:** Safety glasses
- **RESPIRATORY PROTECTIONS:** Work in a fume hood.
- **GENERAL HIGIENIC MEASURES:** The usual precautionary measures when handling chemicals.

9- Physical and chemical properties

- **INFORMATION ABOUT BASIC PHYSICAL AND CHEMICAL PROPERTIES:**
- **APPEARANCE:** Solid
- **COLOR:** White
- **MOLECULAR WEIGHT:** 529,7 g/mol
- **SOLUBILITY:** Good in MeOH, BuOH, fairly in Acet, EtOAc, hardly in Chl, Hex, Water
- **FLASH POINT:** 510.8°C
- **EXPLOSIVE PROPERTIES:** Not determined.

10-Stability and reactivity

- **REACTIVITY:** Not determined
- **CHEMICAL STABILITY:** Not determined
- **THERMAL DECOMPOSITION:** Not determined
- **CONDITIONS TO AVOID:** Heat, flames and sparks.
- **HAZARDOUS DECOMPOSITION PRODUCTS:** Irritating and toxic fumes and gases. Carbon, Nitrogen and Sulphur oxides.
- **INCOMPATIBLE MATERIALS:** Not determined.

11-Toxicological information

- **RTECS#:** ON7050000
- **LOCAL EFFECTS:** Not determined
- **SKIN IRRITATION:** Toxic by skin contact
- **ACUTE TOXICITY:** : DL50 Ip - mouse - 3,1 mg/kg. May affect to liver, kidneys, bladder and intestine.
- **CHRONIC TOXICITY:** Not determined
- **OTHER TOXIC EFFECTS ON HUMANS:** Not determined
- **CARCINOGENIC EFFECTS:** Not determined
- **MUTAGENIC EFFECTS:** Not determined
- **REPRODUCTION TOXICITY:** Not determined
- **TERATOGENIC EFFECTS:** Not determined

12-Ecological information

- **ECOTOXICITY:** No further relevant information available.
- **TOXICITY OF THE PRODUCTS OF:** No further relevant information available.
- **BIODEGRADATION: (Persistence and degradability)** No further relevant information available.

13-Disposal considerations:

- **METHODS OF DISPOSAL: WASTE OF RESIDUES:** Dissolve in or mix the material with a combustible solvent and burn in a chemical incinerator equipped with an afterburner and scrubber.
- **CONTAMINATED PACKAGING:** Disposal must be made according to official regulations.

14-Transport information

- **UN-Number**
- **UN proper shipping name** ADR/RID: Not dangerous goods
IMDG: Not dangerous goods
IATA: Not dangerous goods
- **Transport hazard class(es)** Not applicable
- **Packing group** Not applicable
- **Environmental hazards:** Not determined
- **Marine pollutant:** Not determined

- **Special precautions for user:** Not applicable.
- **Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code:** Not applicable

15-Regulatory information



- **Chemical safety assessment:**

R23/24/25 Toxic by inhalation, ingestion and skin contact.

S36/37/39 Use suitable gloves, eyes and face protection.

S45 If swallowed, seek medical advice immediately. 16

16-Other information

It is the user's responsibility to determine the suitability of this information for adoption of necessary safety precautions. ***This product is intended for research use only*** and should be handled and used only by experienced personnel. Biomar Microbial Technologies shall not be held liable for any damage resulting from the use of this product.