

Parque Tecnológico de León. Parcela M-10.4
24009 Armunia (León) SPAIN
Telephone: +34 987 849 200

León 21/03/2018



Product Code: JU0045
Name: Nocardamine; Norcardamin; Deferrioxamine E; Desferri-ferrioxamin-E; Proferrioxamine-E

1. General information

Molecular Formula:	C27H48N6O9
Molecular Weight:	600,7
CAS #:	26605-16-3
Long-term Storage:	Keep in dark, dry and cool place -20 °C
Assay (HPLC-UV, %)	99,6
Appearance	white powder
Date of Last Analysis	16/02/2017
Retest date	16/02/2022

2. Activity according to literature

Actividad descrita frente a Mycob. Sp.

3. Format

This compound is available in every format from 1 up to 5 mg. In case of needing more than 5 vials handling costs could be applied. Please feel free to send a request for amounts above 5 mg.

For some of our products grams might be immediately available

3. Advice and remarks

Whenever the container is opened for removal of aliquot portions of the substance, the person handling the substance must assure that the integrity of the substance is maintained and proper records of all its handlings are kept. Special care has to be taken to avoid any contamination or adulteration of the substance.

The pharmacological and toxicological properties of this product have not been fully investigated.

Analytical HPLC conditions by request.

Use caution when handling.

For research use only. Not for human or drug use.