

Parque Tecnológico de León. Parcela M-10.4
24009 Armunia (León) SPAIN
Telephone: +34 987 849 200

León 13/12/2018



Product Code: CL0119
Name: Urdamycin A/Kerriamycin B

1. General information

Molecular Formula:	C43H56O17
Molecular Weight:	845
CAS #:	98474-21-6
Long-term Storage:	Keep in dark, dry and cool place -20 °C
Assay (HPLC-UV, %)	96,24
Appearance	Brown amorphous solid
Date of Last Analysis	21/11/2016
Retest date	21/11/2021

2. Activity according to literature

cytotoxic, antibacterial, antibacterial, g-pos., murine I1210 leukemia, ehrlich ascites carcinoma

3. Format

This compound is available in every format from 1 up to 5 mg. In case of needing more than 5 vials handling costs could be applied. Please feel free to send a request for amounts above 5 mg.

For some of our products grams might be immediately available

3. Advice and remarks

Whenever the container is opened for removal of aliquot portions of the substance, the person handling the substance must assure that the integrity of the substance is maintained and proper records of all its handlings are kept. Special care has to be taken to avoid any contamination or adulteration of the substance.

The pharmacological and toxicological properties of this product have not been fully investigated.

Analytical HPLC conditions by request.

Use caution when handling.

For research use only. Not for human or drug use.